



trio ostomy care USA[®]



Trio Ostomy Care USA
Davenport, FL USA

T: 877-400-1700 E: mail@trioostomycare.us



TrioOstomyCare



@TrioOstomyCare



@TrioOstomyCare

trioostomycare.us

trio ostomy care™





Contents

- 4 Silicone is the way of the future
- 7 What is Trio Responsive Silicone®?
- 8 Trio **Siltac**® - Silicone ostomy seal
- 10 Trio **Silvex**® - Silicone convex seal
- 12 Trio **Silex**® - Silicone flange extender
- 14 Trio **Silken**® - Silicone stoma gel
- 16 Trio **Elite**® - Sting free adhesive remover
- 18 Trio **Elisse**® - Sting free skin barrier
- 20 Trio **Pearls**® - Gelling and odour control sachets
- 22 Product codes

All silicone protection... sometimes your skin needs a little help.

Skin is the largest organ of your body and usually provides an excellent barrier against external irritants. For many ostomates however, keeping skin healthy becomes a daily challenge.

Constant contact with adhesive chemicals, coupled with daily removal and exposure to stoma output leaves the skin 'at risk' of damage. Even pink inflamed skin can cause the ostomy flange not to stick, leading to further leaks and a worsening condition.

Trio's unique Responsive Silicone® technology offers many benefits over traditional hydrocolloid ostomy products.

Our unique silicone adhesives:

- Are designed to adhere securely without irritating the skin
- Provide a complete barrier between the skin and the irritant to protect it from damage
- Can be removed clean without the use of adhesive removers but more importantly without stripping skin



trusted results in ostomy

The handy icon guide for our products

Below is a summary of all the icons and their meanings to each Trio product covered in this brochure.



No melting



Flexible



Extended wear time



Protective



Secure adhesion



Sprays at any angle



Self moulding



Waterproof



Odour control



Pain free



No residue



Reduced leakage



Breathable



Silicone technology

What is Trio Responsive Silicone® Technology?

Trio Responsive Silicone® Technology is a unique adhesive designed to protect the skin from harm. It does not absorb the body waste which means it doesn't breakdown or melt, giving your skin the fighting chance it needs to stay safe or recover.



Does not dissolve so provides a durable, long lasting seal
– extends wear time of your ostomy appliance



Reduced risk of leakage by providing a closer, firmer fit
around the stoma



Easily stretches and moulds itself around your stoma,
conforming to your unique stoma shape and size



One piece removal – leaves no residue and no trauma to skin



Provides active protection by responding to body
movement and ensure skin stays safe

Does your **seal** become

Gooney

and **mushy?**



If your flange or seal goes gooey around your stoma, it is because the body waste is being absorbed and melting it. Any skin around the stoma that is exposed will also be 'at risk' of damage.

Common concerns:

- your flange or seal gets gooey and sticky close to the stoma
- you are getting leaks
- your skin is red or sore under your flange
- you have scars or dips around your stoma causing leakage
- your skin is sensitive to your appliance

Main cause for concern:

Your flange and seal are made from a material known as hydrocolloid. This absorbs stool or urine and keeps it against the skin which causes the skin to get red or sore, signs that the skin is becoming damaged.

 **Siltac**[®]



Trio **Siltac**[®] with **Trio Responsive Silicone**[®] is a unique soft silicone ostomy seal that is designed to protect the skin immediately around the stoma from irritation and breakdown by providing a snug fit. It helps to close the gaps around the stoma protecting the skin from stool or urine. Trio **Siltac**[®] doesn't go gooey and mushy like your usual conventional products.

Trio **Siltac**[®] works by responding to the stoma's natural shape, bending, and moving as your stoma works through the day. You simply stretch it around the stoma to provide a secure seal without the need to 'mould' to fit; it just returns to its natural shape making it quick and easy to use.

Key benefits:

- Snug fit around your stoma - eliminating gaps and leaks
- Responds to body movement for a comfortable secure fit
- Doesn't melt or breakdown - providing longer wear time
- Protects skin from stool and urine
- Allows the skin to heal and stay healthy
- Is more hygienic - it does not absorb stoma output



Experience has shown that stomas are as different as the ostomates themselves and they all have their unique needs. With this in mind, Trio has designed Trio Silvex® with Trio Responsive Silicone®, a soft silicone convex seal for use when mild convexity is required.

Over time your body shape will change and so will the requirements of your ostomy management, Trio Silvex® is a convex seal that puts gentle pressure on the skin around the stoma allowing the output to be channelled into the pouch whilst protecting the skin from irritation. The convex shape can fill any dips and 'moats' around the stoma that can cause leaks and sore skin, and it can also help to put gentle pressure on the abdomen to just 'get everything out of the way'.


Trio Silvex® simply stretches around your stoma filling creases and dips whilst providing a firm fit around your stoma reducing any gaps that expose skin which is at risk of breakdown and irritation. Trio Silvex® removes in one piece leaving no residue – your pouch change regime can now be quicker, easier and cleaner.

Key benefits:

- Provides gentle pressure to form a snug fit around your stoma
- Responds to body movement
- Doesn't melt or breakdown - providing longer wear time
- Protects skin from stool and urine
- Allows the skin to heal and stay healthy


Hydrocolloid

Gaps around the stoma expose delicate and sensitive skin.



Trio Responsive Silicone®

Returns to its natural shape around the stoma, which helps to remove gaps and protect your skin.



Does your **convex seal**
give you extra

SUPPORT?



Does the
edge of your bag

Lift?





Ostomates often need additional security around their flange for a variety of reasons:

- Peristomal hernia
- Physical activity
- Body shape
- More surface area
- Increased security especially during activity and at night

These are just a few common reasons that Trio **Silex**® soft silicone flange extender is fast becoming the product of choice. Ultra-thin **Silex** is designed to go around the flange to provide a soft, flexible border which conforms to your body shape to provide a security even during sleep.



Trio **Silex**® is waterproof which is ideal for bathing, swimming, and warm conditions as it does not absorb moisture. This means it will **not breakdown or 'melt'** and can be **gently peeled away** from the skin without the need of an adhesive remover, leaving residue or damaging the skin.

Key benefits:

- Secure adhesion during activities, warm conditions and varied body shapes
- Ultra thin, invisible skin for comfortable instant wear
- Waterproof – ideal for swimming, bathing and holidays
- Highly flexible – allows for more freedom of body movement
- Gentle removal without any pain and leaves no residue
- Allows the skin to breathe – reduced risk of maceration and healthier skin



There are times that your stoma just doesn't want to play ball, or your skin is sore or has creases that causes leaks. In such situations, you may require Trio Silken® silicone ostomy gel to seal the tiniest gaps that allow output under the flange and cause leaks and sore skin.

Why is Silken different to other pastes?

- It is made from silicone which does not 'melt' into a gooey mess
- It does not contain alcohol or preservatives so it won't sting or irritate the skin
- It does not require air or heat to set – it uses the skin's own moisture to turn into a 'rubbery' seal
- Trio Silken® works best when allowed to settle into any dips and creases as you move around after applying your pouch. It will collect in these dips and creases, and set to close any channels that allow leaks
- Less is more! Only a small amount is needed where the leak usually occurs. No need to spread, simply apply your pouch in the normal way.
- It comes off clean from the skin during your routine bag change without the need for adhesive removers



Trio Silken® is a silicone gel that is used to fill skin folds, dips and creases, and uneven areas around the peristomal skin that causes the output to seep under the flange. Maintaining healthy skin around your stoma is the secret to good stoma management and Trio Silken® is the key.

Does your
paste

Stretch

with you?



Do you experience

painful

removals?





The next generation of silicone medical adhesive removers has arrived! A new 100% silicone formulation that has been carefully designed for you. Trio Elite® does not use a traditional propellant and will not feel cold when coming in contact with your skin. Spraying at any angle, with a gentle action, Trio Elite® has a controlled & directional spray pattern for longer lasting use Elite by name, Elite by nature!

Key benefits:

- Will not sting, even on sore skin
- Releases adhesives quickly
- Sprays at any angle
- Controlled & directional spray pattern for longer lasting use
- Will not affect adhesion of the next Ostomy appliance, dressing or tape
- Dries in seconds – leaving no residue
- Does not look like a medical adhesive remover

Wound Care - dressings and fixation devices

The new, 100% silicone technology of Trio Elite® will painlessly remove dressings, tapes and fixation devices without causing trauma to the skin.



Ostomy Care - pouches

Ensuring a good seal between the ostomy wafer and the skin is essential to preventing leakage. The strong adhesives necessary to achieve this can cause pain and trauma on removal often damaging your peristomal skin. Trio Elite® can release even high-tack adhesives quickly & efficiently!

Remove unwanted sticky adhesive residue

Most medical adhesives are designed to adhere to undamaged, clean skin. Gently remove unwanted sticky residues with Trio Elite® soft wipes to ensure the skin is prepared for the next ostomy appliance or adhesive product.



trio Elisse®



Protecting the skin from damage and having healthy peristomal skin is the key to good ostomy management. We are not designed for things to stick to us; we sweat, move, produce heat and oils, but if the skin is damaged then this adds a new level of complication. Trio Elisse® silicone skin barrier acts as the perfect defence against attack from body waste and chemicals from adhesives, whilst allowing the skin to stay safe underneath.



Trio Elisse® is the **world's first ever skin barrier** to have a patent due to the uniqueness of the formulation.

Trio Elisse® creates a breathable and durable barrier between the skin and the wafer and its hypoallergenic formulation soothes sore and reddened skin and actively reduces itching.

Key benefits:

- Protects skin from damaging effects of adhesives and body fluids
- Unique patented formulation soothes reddened and sore skin
- Will not sting even on sore skin
- Dries quickly
- Sprays at any angle
- Does not cause build-up of layers
- Manages moisture on mildly damaged skin



Does your pouch
make your skin

rawl?



Do you **dream** of a

Peaceful

night's sleep?





Trio Pearls® Super absorbent gelling and odour control sachets have been carefully formulated to provide ostomates who have loose stool, gain more control over the management of their ostomy and their lives.

Trio **Pearls®** can provide control by absorbing the fluid in loose stool, and gelling inside the pouch. The formed gel is thicker in consistency which helps to prevent leaks normally associated with liquid output; prevents the filter from clogging prolonging its efficiency; and avoids bulking by allowing the content to be flattened.

Active Odour Control

Trio **Pearls®** is not just a gelling agent – it has an organic compound which specifically targets odour.

The breakdown of certain foods in the digestion process produces odorous molecules and in turn gases. Ostomates often fear the loss of control over their bodily function following stoma surgery, and flatus and odour can be a major concern affecting quality of life for ostomates.

NanoActive Technology

The new Trio **Pearls®** incorporate a superabsorbent gelling crystals and nano-engineered ingredients that actively adsorb odorous organic molecules rather than just trapping them in the gel.

Key Benefits:

- Actively adsorbs odour
- Preserves filter
- Prevents leakage
- Make the pouch more discreet
- Easier drainage/less messy
- Gel does not stain clothes or toilet
- Security over longer periods

Perfect for:

- Overnight use for a good night's sleep
- Long journeys or flights – take it on holiday
- Social occasions such as parties, restaurants, theatre/cinema, etc.
- Active pursuits
- Everyday use

Available now - Ordering information

Code	Description	Pack	HCPCs	Pages
TR1020	Trio Siltac® 1 - Silicone Ostomy Seals (20mm-28mm)	30	A4385	8-9
TR1028	Trio Siltac® 2 - Silicone Ostomy Seals (28mm-35mm)	30	A4385	8-9
TR1035	Trio Siltac® 3 - Silicone Ostomy Seals (35mm-44mm)	30	A4385	8-9
TR1020CX	Trio Silvex® 1 - Silicone Convex Seals (20mm-30mm)	10	A4411	10-11
TR1030CX	Trio Silvex® 2 - Silicone Convex Seals (30mm-40mm)	10	A4411	10-11
TR1060	Trio Silex® Silicone Flange Extenders	20	A4362	12-13
TR1070	Trio Silken® Stoma Gel	60g tube	A4406	14-15



Code	Description	Pack	HCPCs	Pages
TR201	Trio Elite® Sting Free Adhesive Remover Spray	50ml	A4455	16-17
TR202	Trio Elite® Sting Free Adhesive Remover Wipes	30	A4456	16-17
TR203	Trio Elisse® Sting Free Skin Barrier Wipes	30	A5120	18-19
TR204	Trio Elisse® Sting Free Skin Barrier Spray	50ml	A4369	18-19
TR205	Trio Pearls® Gelling & Odour Control Sachets	100	A4422	20-21



#ONLYWITHTRIO

Trio Healthcare are the only company that have brought ostomates a new solution for skin problems and leaks in almost 40 years.

Trio® Healthcare is committed to making a positive difference in the care of ostomates by making their products cost effective and efficient. If you like Trio® products and wish to have them on your insurance, you can discuss it with your retailer. If you have any questions about Trio® products, give us a call and we would be pleased to help you.

